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Hand Stretching Rehabilitation Therapy Device

CATALOG

Ziglite Smart Health Care Products Co., Ltd.
節亮康護產品有限公司
Your Smart Health Care Solution Partners

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Hong Kong Office Address

Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road, Tsuen Wan, Hong Kong

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
Ziglite Smart Health Care Products Co., Ltd. Incorporated in 2018, we work with our partners to bring their products to the local market. Mainly engaged in the management of health care products, the main customers are hospitals and clinics, homes for the elderly and home patients.

We are committed to the development of excellent health care products, dedicated to the needs of a wide range of different customers, get the best care in the family treatment process, and greatly improve the effect of rehabilitation. Get the best care during the family treatment process and greatly improve the healing effect. We will provide quality products with stable quality and stable performance for customers.

We specialize in a wide range of products including: anti-virus and air freshening, incontinence care, mobility aids, decompression supplies, etc.

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Hand Stretching Rehabilitation Therapy Device

手指選擇
Select finger

循環訓練
Cycle training

鏡像訓練
Mirror training

ZIGLITE[®]
HEALTHCARE
SBA-2

循環訓練+鏡像訓練+手指選擇訓練
Cycle Training + Mirror Training + finger Selection Training

- 手部伸展復康治療器在患者整個恢復週期均可參與，通過對患者患肢的訓練，刺激患者神經、肌肉，協助患肢恢復功能，提高患者生活質量。
- 適用於受傷後患者的手指關節恢復或由腦損傷、腦卒中、腦癱、脊髓損傷、骨科術後、多發性硬化等引起的肢體障礙後遺症患者在家庭及機構中的康復訓練。
- The hand stretching rehabilitation device can participate in the patient's entire recovery cycle. Through training the patient's affected limb, it stimulates the patient's nerves and muscles, assists the affected limb to restore function, and improves the patient's quality of life.
- It is suitable for the rehabilitation of the finger joints of patients after injury or the rehabilitation training of patients with sequelae of limb disorders caused by brain injury, stroke, cerebral palsy, spinal cord injury, orthopedic surgery, multiple sclerosis, etc. at home and in institutions.

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Product Introduction

Training Terminal

- Principles of human kinematics
- Flexible flexion and extension
- Stretching / fist clenching: large tilt angle
- Enhance joint mobility and muscle strength recovery, and promote active flexion and extension movement

Rehabilitation Therapeutic Device Host

- cycle training, mirror training are optional and finger selection training
- 5 levels of power adjustable
- The default is a 30-minute training countdown, and a maximum of 60 minutes of training can be set


Mirror Training Gloves

- Mirror linkage training
- The healthy hand drives the affected hand
- Real-time pressure sensing

Application Scenario

Rehabilitation institutions

Pension institutions

Patient home training

Product Features

Simulate Muscle Drive

The driving module is fitted with the human joint position, and the degree of freedom of movement is the same as the degree of freedom of hand joint movement. It simulates the human kinematics mechanism. The driving module can be dynamically adjusted according to the patient's movement state to avoid secondary injury to the patient.

Circulation Training To Relieve Muscle Tension

The circulation training mode adopts fully automatic design, which helps to promote the recovery of hand muscle strength, increase joint range of motion and promote the emergence of active flexion and extension movement.

Mirror Training Enhances Control

Mirror training can map the movement pattern of the healthy hand to the affected hand, and stimulate mirror neurons through mirror movement to promote the remodeling and recovery of brain function and hand function.

Finger Selection Training

Finger selection training can carry out rehabilitation treatment for single or multiple fingers according to the actual needs of patients, which is more targeted.

Technology Driven

Mirror training and passive training are combined to promote rehabilitation

Finger Selection Training

Cycle Training

Mirror Training

LCD Display

According to the principle of brain reconstruction, task-oriented training is combined with hand rehabilitation robot training, Help patients complete controlled motor function training, Through repeated intensive training, it can promote the functional reconstruction of the brain central nervous system.

Product Introduction

Training Gloves

- Using the principle of software robot technology, flexible drive, avoid causing secondary injuries such as strains to the patient.

Rehabilitation Treatment Device Host

- 5 levels of strength selection, to help you quickly find the training intensity you need.
- The smart button can set the training time period at will.

Cycle Training


- Fully automatic design is adopted to promote the flexion and extension of affected fingers, muscle strength recovered and joint range of motion increased.

Mirror Training


- The healthy hand moves the affected hand, visual stimulation activates mirror neurons, promotes the emergence of brain active motor function, and solves the basic rehabilitation problems of hand patients in the early stage of rehabilitation period at will.


Training Gloves



Mirror External Control Gloves

Power Adapter

Training Host

Function list of therapeutic device operation buttons

Number	Button	Function Description
1		Press the "Cycle Mode" button to switch to the cycle mode, in this mode, press again to enter the finger movement training selection interface. Long press the "Cycle Mode" button for 3 seconds and release it to set wave finger training and all mode selections. In the finger movement training selection mode, you can select / pause thumb training..
2		Press the "Mirror Mode" button to switch the mirror mode. In the finger movement training selection mode, you can select / pause the index finger training.
3		In finger movement training selection mode, click "middle finger cycle selection" to select / pause middle finger training.
4		Each time you press the "Training Strength" button once, you can set the strength of the training. The set values are: 1, 2, 3, 4, and 5, with a total of 5 strengths. In finger movement training mode, ring finger training can be selected / suspended.
5		Each time you press the "Training Cycle" button once, you can set a grip and stretch training cycle time. The set values are: 1, 2, 3, 4, and 5, with a total of 5 settings. In finger movement training mode, you can select / pause ring finger training.
6		Each time you press the button once, you can start, pause, and return in the work training mode.

Special function settings

Number	Function	Functional operation description
1	Circuit training Time setting	In the main interface of the cycling training mode, long press the "⌚" button for 3 seconds to enter the time setting mode, and the exercise time digits will flash. The "📊" and "⌚" button can be used to set the addition and subtraction of the training time. The minute setting division value is 5 minutes, and the training time can be up to 60 minutes.
2	Mirror mode Strength setting	In the mirror image mode, each press of the "📊" button can set the addition or subtraction of the training intensity. There are 5 settings in total.
3	Mirror mode Shrink settings	In the mirror mode, press the "⌚" button once each time to set the length of the training contraction cycle, a total of 5 gears.

Buttons And Screen Icons In De

Operation buttons and interface display

Therapeutic device operation interface display list

Number	Icon	Function description
1	⌚30:00	The training time of the therapeutic device is 30:00 minutes by default and can be set to 60:00 minutes at most. The training begins to count down, and the therapy device automatically stops after returning to zero.
2	力量:📊	It is a display of the training strength of the treatment device, and the strength can be set to 1, 2, 3, 4, and 5 levels of strength.
3	週期:📊	The training cycle of the therapeutic device is displayed, and the cycle can be set to 1, 2, 3, 4 and 5.
4	收縮:📊	The image training contraction display of the therapeutic device shows that the operation time of the contraction motor can be set as: 1, 2, 3, 4 and 5.
5	動作👉	Selection and setting of finger exercises. Settable range: thumb, index finger, middle finger, ring finger, and tail finger.

Detailed button function

The diagram shows the device interface with numbered callouts explaining button functions:

- 1) Cycle / split finger mode switching Thumb training options:**
 - Press the button once to enter the finger selection interface: corresponding to the finger button to select/stop a finger training.
 - In the finger selection interface, press and hold the button for 3s to enter the action mode switch: ①Wave cycle training for each finger/②Cycle training mode switch for all fingers.
 - In the main interface, press and hold the button for 3 seconds to enter the training time setting, up to 60 minutes.
- 2) Mirror mode Index finger training:**
 - Press the button once to switch the mirroring mode.
 - In the finger selection interface, click to select/stop the index finger training exercise.
- 3) Middle finger training options:**

In the finger selection interface, click "Middle finger cycle selection" to select/stop the middle finger training exercise.
- 4) Power setting Ring finger training options:**
 - Press the "Training Strength" button to set the strength of the training, the set values are: 1, 2, 3, 4, 5, a total of 5 levels of strength.
 - In the finger selection interface, click Select/Stop ring finger training.
- 5) Cycle time setting Tail finger training options:**
 - Press this button to set the cycle time of one grasping and stretching training. The setting value is 1 ~ 5, a total of 5 gears.
 - In the finger selection interface, click optional / stop tail finger training.
- 6) Start/pause Return:**

Start, pause, return and other operations can be performed in the work training mode.


» Interface display details

Cycle mode display

Other interface display

● Mirror mode display


● Select finger display

(Flashing means selected, stop flashing means to pause training)

1 Each finger wave type circuit training mode

2 All finger loop training modes

User Guide

How to use the hand stretch rehabilitation therapy device?

1. Cycle training

The circulation training mode adopts a fully automatic design. When the patient has no active control ability, it is recommended to use circulation training first, which will help to promote the recovery of flexor and extensor muscle strength of the affected hand, increase the joint's range of motion, and promote the emergence of active finger extension movement.


» Default circuit training

<p>1</p>  	<p>1) Wear Training Gloves 2) After startup, enter the main training interface</p>	<p>2</p>  <p>3) After powering on, press the start button to enter the default training: split-finger wave cycle mode, training time is 30 minutes, cycle strength is the third gear.</p>
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

» Set up circuit training as needed

1 Power, period, time setting (4 and 5 can be skipped if you don't need to set the time)

<p>(1)</p> 	<p>(2)</p> 	<p>(3)</p> 	<p>(4)</p> 
<p>Wear training gloves Enter the main training interface after booting.</p>	<p>Set training intensity Press 1 time to set the gear, 1 ~ 5 gears in total.</p>	<p>Set training cycle Press 1 time to set the gear, 1 ~ 5 gears in total.</p>	<p>Set training duration Press and hold for 3 sec to enter the setting, and set time plus and minus, each click changes 5 min.</p>

(5) After setting, click to return to the main interface.

2 Training finger selection (Finger icon flashes for training state, static state is for stop state)


<p>(6)</p> 	
<p>Training finger selection In the main training interface, click to enter the finger selection interface, select training/stop a finger training.</p>	<p>Finger control: select/stop switching each time you press In this mode, the buttons from left to right correspond to thumb, index finger, middle finger, ring finger and tail finger control.</p>

3 Training loop mode selection

<p>(7)</p>  <p>動作 Finger wave loop 動作 Finger circulation</p>	<p>(8)</p>  <p>Return Start training</p>
<p>Training finger selection In the main training interface, long press for 3 sec to enter the action mode switch: ① Wave type circuit training for each finger ② Common circulation training mode for all fingers</p>	<p>Start circuit training After setting, click to return to the main interface. Click again to start training.</p>

2. Mirror Training

Mirror mode refers to the mode that the patient drives the patient's hand to move synchronously through healthy hand movement. Through interactive visual feedback, it can stimulate different regions of the patient's brain to produce motor and sensory consciousness, start mirror neurons, and promote the reconstruction of motor function and brain function.

<p>(1)</p> <p>Training glove Mirror glove</p>	<p>(2)</p>	<p>(3)</p>	<p>(4)</p>
<p>Put on gloves and turn on</p> <p>Wear training gloves on the affected hands, Wear mirror gloves for healthy hands , After startup, enter the main training interface.</p>	<p>Enter mirror mode</p> <p>In the main training interface, press the button once to enter the mirror training mode.</p>	<p>Set training intensity Set shrink duration</p> <p>Press once each time to set, and the training intensity is 1 ~ 5 gears in total.</p> <p>Press once each time to set, and the contraction duration is 1 ~ 5 gears in total.</p>	<p>Start mirror training</p> <p>The healthy hand grasps (extends) the movement to issue the instruction, and the affected hand grasps (extends) the movement according to the instruction.</p> <p>Note: you need to grasp and press (release) the control button in the mirror glove to start the grip (release) action of the training glove.</p>
<p>Control the corresponding finger training:</p> <p>There are 5 control keys in the mirror glove, and each key corresponds to 5 fingers of the training glove. The healthy hand sends movement instructions to the corresponding fingers of the training gloves by controlling the buttons in the mirror gloves, and the corresponding fingers of the training gloves receive instructions to make grasping (extension) movement</p>			

Technical Parameters

Product specification

Product Name	Hand Stretching Rehabilitation Therapy Device	Liquid inlet protection degree	IPX0
Product Model	SBA-2	Strength setting	1-5 files
Security Type	Type I	Time setting	5-60 minutes
Power Supply	AC100-240V, 50/60Hz	Operating noise	≤60dB
Host Size	251*204*111mm	Training mode	Cycle training, Mirror training


Glove size chart

Size	S	M	L	XL	XXL
Hand length (mm)	140-160	160-180	165-190	190-220	220-230
Hand width (mm)	60-80	70-90	80-100	100-120	120-130
Middle finger circumference (mm)	45-50	50-60	60-70	70-90	90-100

* The corresponding model can be selected according to the user's actual hand measurement situation.

If two or more measured values are at the critical point, it is recommended to choose one yard larger.

* Women generally choose s/m, men generally choose l/xl, and the larger hand is XXL.

! Safety Precautions

[Warning]

Failure to follow the following safety instructions may lead to serious accidents such as fire, electric shock and casualties!

- Do not disassemble, repair or modify the components of the equipment by yourself to avoid causing abnormalities of the therapeutic device.
- Do not store components inside or near the water heater, microwave oven, cooking equipment or high-pressure container to avoid fire caused by overheating.
- Do not sterilize or autoclave.
- Do not continue to use when there is an alarm prompt.
- It is suitable for hand circulation and mirror motion rehabilitation training. Do not use it for other purposes to avoid accidents.
- Do not store all parts in an overheated or supercooled place to avoid accidents.
- Please keep the power connector dry without water. Please do not plug and unplug the power with wet hands to avoid electric shock or injury.

I Storage and transportation

- During transportation, please fix the therapeutic device in a relatively stable place to prevent collision and damage of the therapeutic device.
- Please prevent liquid or foreign matters from entering the therapeutic device and avoid damaging internal electronic components and causing faults.
- Do not store in the environment of disinfectant gas or active gas to avoid damaging internal electronic components and causing faults.
- Do not place it in strong light irradiation places such as sunlight or ultraviolet for a long time to avoid discoloration of the shell.
- Do not place this product in a place accessible to infants and young children.
- When saving the product, do not place other items on the product.
- Do not let the product touch nails, needles and other sharp objects.
- Do not fold the gloves forcefully to avoid damaging the software actuator and causing failure.
- Do not store the therapeutic device in a humid environment or strong magnetic field interference, which is not suitable for the working environment of electronic equipment.

I Working environment

- Ambient temperature: 10 °C ~ 35% °C
- Relative humidity: 30% ~ 75%, excluding condensation
- Atmospheric pressure range: 700hPa ~ 1060hpa
- Power supply: AC100-240V, 50 / 60Hz (therapeutic device host power supply DC12V)

I Post processing

- Do not destroy the therapeutic device in the fire.
- Do not treat the therapeutic device as ordinary domestic waste. At the end, it shall be treated according to the requirements of medical device standards. Do not discard it at will.

Safety Precaution Include

- In case of skin allergy and other discomfort during the use of the product, please stop using and consult a doctor.
- For patients with diseases or wounds or being treated for them, or with previous medical history, please consult a doctor after being familiar with the following contents.
- If there is sprain, trauma, bone injury or fracture of unknown condition on the hand, people who have not been treated or need braking after treatment should not use this product without professional confirmation.

⊞ Contraindication

- Patients with severe blood flow disorders such as thoracic venous aneurysm and hemocolumn disease.
- Patients with ulcer, wound, purulent fruit, bleeding, scald, inflammation, etc. on the affected hand.

⊞ Other Precautions

- If you experience an abnormality during use or your body feels abnormal, you should stop using it and use it under the guidance of a doctor after consulting a doctor.
- Redness, itching, rash, etc. appear during use, stop use and seek medical attention.
- Do not use with other equipment at the same time or use it while applying medicine.
- Do not use it alone by unconscious people or young children.
- Do not use this product for purposes other than the purpose of use.

Before Use

- 1) Before wearing training gloves, please try to run the therapeutic device host for 1 minute to confirm that all parts function normally.
- 2) For patients with mandatory systolic spasm, massage before using the therapeutic device.
- 3) Note: before training, please confirm that the non muscle tension of the hand to be trained is too high (above grade II).

In Use

If the following conditions occur during operation, please stop using it immediately and disconnect the power supply!

- 1) Abnormal sound from the host
- 2) The surface temperature of the host is too high
- 3) The host accidentally enters water
- 4) The user feels pain, swelling or other discomfort
- 5) If spasm occurs during training, stop using it immediately

Usage Time

- 1) It is suggested to select appropriate training methods according to their own actual situation
- 2) The time of each use should not be too long. If the hand feels sore the next day after training, it is necessary to reduce the number and time of daily training
- 3) It is recommended to use 2-4 times a day for about 20 minutes each time, with an interval of at least 45 minutes
- 4) If you need to change or extend the use time, please consult professional medical personnel

⊞ Use Of Power Adapter

- The power adapter must be plugged in and out by the patient's supervisor.
- After using the therapeutic device, please cut off the power adapter and 220V power supply.
- Do not use a power adapter that is not provided with this machine.