

Lifting Ring Stretching For Upper Limb Rehabilitation Training

上肢康復訓練吊環拉伸

Model : RP-IX

Features:

- Pulley alloy steel material is more durable, sponge handle, wear-resistant tension rope, comfortable and durable
 - Upper limb rehabilitation training, 360 ° Rotation
 - Applicable population:
 - 1) All patients with upper limb dysfunction (including shoulder joint dysfunction) (including stroke, cerebral infarction, cerebral hemorrhage, hemiplegia, etc.)
 - 2) Shoulder and neck disease patients, elderly rehabilitation exercise, white-collar workers
-
- 滑輪合金鋼材質更耐用, 海綿手柄, 耐磨拉力繩舒適耐用
 - 上肢康復訓練, 可360°旋轉
 - 適用人群:
 - 1) 一切存在上肢功能障礙 (含肩關節障礙) 的病人 (包含中風、腦梗、腦溢血偏癱等引起的)
 - 2) 肩周炎、頸椎增生及病變的病人, 老年人康復健身鍛煉使用, 長期伏案的白領工作人群



Specifications

| | |
|-----------------------|---|
| Product Model 產品型號 | RP-IX |
| Product Function 產品功能 | upper limb exercise 上肢鍛煉 |
| Product Material 產品材質 | Aluminum alloy pulley + sponge handle + wear resistant tension rope 鋁合金滑輪+海綿手柄+耐磨拉力繩 |
| Product Size 產品尺寸 | 30*4cm |
| Product Weight 產品重量 | 238g |
| Packing Size 包裝尺寸 | 13*12.5*8.5cm |
| Package Weight 包裝重量 | 307g |

Instructions

- ◆ 首先將背帶調整到合適的長度
 - ◆ 穿上後打開開關, 然後站立或者坐正
 - ◆ 保持4秒後, 感應器會震動一下, 感測器將您的正確姿勢記錄下來
 - ◆ 當您背部有彎曲弓腰駝背時, 會震動提醒您挺胸抬頭
-
- ◆ 設置:
長按3秒震動一次, 藍色燈亮, 表示機器已開機
姿勢在4秒內保持不變, 此時連續震動兩次表示角度已設定, 校正角度為設定角度+偏移量25度報警
再次長按3秒, 震動一次, 藍色燈熄滅, 表示關機